

**No water, no mess.  
Unlike the rest.**

Thanks to its unique, self-hydrating coating, MAGIC<sup>3</sup> Go<sup>®</sup> with CLEANGLIDE<sup>™</sup> technology provides you with control, comfort and confidence. You never have to deal with the mess of water or lubrication because the coating is already hydrated. Simply open the package, and it's ready to go—no matter where life takes you.



No one chooses to use catheters everyday, but for people who need to, we promise to provide great products, services and support from your doctor's office to your home, your life and anywhere you love to be.

**No water,  
no mess.**



Visit us at [bardcare.com](http://bardcare.com) to join a community for inspiration, resources and support.



This general information is not intended to be a substitute for the advice provided by a qualified health care professional.

*Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.*

8195 Industrial Blvd., Covington, GA  
800-243-3315 | [bardcare.com](http://bardcare.com)

Bard, Magic<sup>3</sup>, Magic<sup>3</sup> Go, and CleanGlide are trademarks and/or registered trademarks of C. R. Bard, Inc.  
©2018 C. R. Bard, Inc. All Rights Reserved.  
1801-32 BMD/INCA/0218/0208



**Magic<sup>3</sup> Go<sup>®</sup>**  
hydrophilic intermittent catheter  
with CLEANGLIDE<sup>™</sup> technology

16" male length catheters



Simply unlike any other catheter.

## Control

## Comfort

## Confidence


### SUREGRIP™ insertion sleeve

The MAGIC<sup>3</sup> Go<sup>®</sup> provides control with a no-touch insertion technique—allowing the catheter to be held firmly without touching the surface, minimizing the risk of contamination.




### New, unique CLEANGLIDE™ technology

CLEANGLIDE™ is a lubricious hydrogel that maintains the catheter in a hydrated state without requiring activation or lubrication prior to use.



### Three-layer silicone catheter

With its self-hydrating coating, the MAGIC<sup>3</sup> Go<sup>®</sup> with CLEANGLIDE™ technology is designed to glide easily and comfortably through the urinary tract.



### Discreet packaging

The MAGIC<sup>3</sup> Go's package is discreet and gives you the ability to carry it folded in a compact format. The adhesive folding tab frees up your hands by attaching to dry, vertical surfaces.

Easy to carry



Easy to open



Easy to use



### Saves time

Thanks to the MAGIC<sup>3</sup> Go<sup>®</sup> catheter there is no need to activate the coating or add lubricant anymore. It is ready to go out of the package.

### Tapered tip and comfort sized eyelets

Enables ease of catheter movement in sensitive areas of the anatomy.

Straight tip

Coudé tip



Request a  
**FREE**  
Sample

Request a free sample and review our available French sizes at [BardCare.com](http://BardCare.com), or call us directly at 800.243.3315.



Not made with natural rubber latex. Does not contain DEHP or phthalate